


# ===== WPI =====

- TI - Bunion correction jig - has press-fit member attached to belt made of cloth-like material
- AB - J10234759 The jig has a press-fit member (1) which is attached with a belt (2) made from a cloth-like material. Support tapes (a,b) are fixed at the belt end parts.
- The press-fit member is inserted between a hallux (5) and a second toe (6) of a human foot (4). The belt is rolled around the hallux and the second toe, and is fixed using the support tapes.
- ADVANTAGE - Gives stabilized force to hallux of bunion. Enables simple detachment. Improves correction effect of bunion.
- (Dwg.2/4)
- PN - JP10234759 A 980908 DW9846 A61F5/10 003pp
- PR - JP970076426 970220
- PA - (SUGA-I) SUGAWARA M
- DC - P32
- IC - A61F5/01 ;A61F5/10
- AN - 98-535384 [46]

# ===== PAJ =====

- TI - VALGUS HALLUX CORRECTING JIG
- AB - PROBLEM TO BE SOLVED: To provide stable corrective force to the hallux of valgus hallux.
- SOLUTION: The device is composed of an insert 1 to be held between the big toe and the second toe, and a belt or cylinder of cloth or expansion material for stably holding it. For adjusting thickness of the insert, a thickness adjusting member is used. In the case of the belt, it is stopped by a surface tape or the like. Simultaneously with the execution of above, a supporter to tighten a middle foot ligament part may be used for further improving correction effect. By thus stably providing corrective force to the big toe of hallux valgus, hallux valgus correcting effect can be improved.
- PN - JP10234759 A 980908
- PD - 98-09-08
- ABD - 981231
- ABV - 098014
- AP - JP970076426 970220
- PA - SUGAWARA MASATOSHI
- IN - SUGAWARA MASATOSHI
- I - A61F5/10; A61F5/01



<First Page Image>